



Custom
Printed
Products



TPJDS100GB

Post-it® DS100 Pop-Up Note and Flag Dispenser Executive Gift Set

8 days production time

Available undecorated or with a custom imprint on the dispenser. Dispenser features a large area for a custom laser engraving or ink imprint.



PDWIGB10

Post-it® Note and Flag Pen Executive Gift Set

7-12 days production time

Available undecorated or with a custom imprint on the pen or pads.

Gift Sets

Post-it® DS330 Pop-Up Note Desktop Dispenser Executive Gift Set

8 days production time

Available undecorated or with a custom imprint on the dispenser. Dispenser features a large area for a custom laser engraving, paper insert or ink imprint.

TPJDS330GB



Post-it® Flag Pen and Flag Highlighter Executive Gift Set

5 days production time

Mix and match any flag pen and flag highlighter for a unique gift set!



WIPHGB

Also Available—
Gift Sets in Support
of Breast Cancer
Awareness!
See page 38.



Post-it® Flag Pen and Flag Highlighter Executive Gift Set										Quantity Price Per Gift Set		
Item #	Description and Imprint	100-249	250-499	500-999	1000-2499	2500-4999	5000-9999		Set-Up	Units/Case Wgt	sm@rt code	
	Includes: Any one (1) Post-it® Flag Pen (WIPFC) and any one (1) Post-it® Flag Highlighter (WIHFC) packaged together in a Silver Gift Box. Each Pen/Highlighter holds 50 flags. Imprint your logo or message on barrel in up to 4 colors.									250 sets/ 25 lbs.		
	1 Color Spot Imprint on Barrel of Both Items	7.20	6.05	5.38	4.90	4.60	4.05	BBCCD	FREE			
	2 Color Spot Imprint on Barrel of Both Items	8.30	7.03	6.26	5.66	5.30	4.27	BBCCD	FREE			
	3 & 4 Color Spot Imprint on Barrel of Both Items	quote	quote	7.46	6.76	6.03	4.59	BCCD	FREE			

Allow 5 days production time. Flags cannot be customized. Flag refills available only with accompanying order, call for quote. Ink and pigment are not refillable. Flag, Pen and Highlighter Combo (WIPHFC) not available for gift set. 4 color pricing does not include 4-color process printing.

Send us your business card or personalized greeting card to be inserted at no extra charge!
See below for availability.

Post-it® DS100 Pop-Up Note and Flag Dispenser Executive Gift Set										Quantity Price Per Gift Set		
Item #	Description and Imprint	10-24	25-49	50-99	100-249	250-499	500-999	1000-1499		Set-Up	Units/Case Wgt.	sm@rt code
	Includes: One (1) DS100 Pop-Up Note and Flag Dispenser, one (1) 3" x 3" 50-sheet Post-it® Pop-Up Note Pad in Canary Yellow (unprinted), one (1) 1" Red flags and one (1) 1" Blue flags packaged together in a Silver and Black Executive Gift Box.										10 sets/ 17.44 lbs.	
	No imprint	17.94	17.94	17.94	16.89	16.89	15.84	15.84	C	FREE		
	1 Color Ink Imprint on Dispenser Top (one or both sides of note opening)	N/A	N/A	23.08	20.86	18.46	16.88	16.67	CCDD	FREE		
	2 Color Ink Imprint on Dispenser Top (one or both sides of note opening)	N/A	N/A	26.76	24.01	20.98	18.87	18.56	CCDD	FREE		
	1 Laser Engraving* (one side of note opening)	N/A	N/A	23.43	21.93	20.13	18.56	16.96	CCDD	FREE		
	Supplied Business/Greeting Card Insert**	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE		

Allow 8 days production time. Actual Post-it® Note size: 2-7/8" x 2-7/8". Unprinted, 3" x 3" x 50- or 100-sheet Post-it® Pop-Up Note Pad refills and Post-it® Flag refills are available retail. *Laser engraving is on underside of dispenser top on one side of opening. For laser engraving in more than one location, call for quote. For both laser engraving and ink imprint on same dispenser, call for quote and production time. ****Supplied insert cards must be sent within three (3) days of receipt of order.**

Post-it® Note and Flag Pen Executive Gift Set										Quantity Price Per Gift Set			
Item #	Description and Imprint	25-49	50-99	100-249	250-499	500-999	1000-2499	2500-4999	5000-9999		Set-Up	Units/Case Wgt.	sm@rt code
	Includes: Ten (10) 3" x 4" x 40 sheet Post-it® Executive Notes with Executive Stock Granite Design and any one (1) Post-it® Flag Pen (WIHFC) packaged together in a Silver and Black Executive Gift Box.											10 sets/ 8.17 lbs.	
	No imprint	6.99	6.99	6.99	6.99	6.99	6.99	6.99	6.99	B	FREE		
	1 Color Spot Imprint on Pen Barrel** add	N/A	1.95	1.15	.70	.50	.40	.35	.30	G	FREE		
	2 Color Spot Imprint on Pen Barrel** add	N/A	2.45	1.30	.75	.55	.45	.40	.35	G	FREE		
	3 & 4 Color Spot Imprint on Pen Barrel** add	N/A	quote	2.30	1.25	.85	.60	.50	.40	G	FREE		
	1 Color Spot Imprint on Notes* add	N/A	2.50	1.25	.55	.25	.15	.10	.08	G	FREE		
	2 Color Spot Imprint on Notes* add	N/A	4.00	2.00	1.00	.40	.30	.18	.12	G	FREE		
Supplied Business/Greeting Card Insert***	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE			

Allow 7 days production time. Custom imprint on notes does not include Executive Stock Granite Design. Flag, Pen and Highlighter Combo (WIPHFC) and Flag Permanent Marker (WIMFC) not available for gift set. *Add five (5) working days production time for custom printed pads. For Executive Stock Granite Design, add one additional imprint color. **Choose any flag pen (WIPFC). Imprint on barrel only. No imprint on flags. 4 color pricing does not include 4-color process printing. *****Supplied insert cards must be sent within three (3) days of receipt of order.**

Post-it® DS330 Pop-Up Note Desktop Dispenser Executive Gift Set										Quantity Price Per Gift Set		
Item #	Description and Imprint	25-49	50-99	100-249	250-499	500-999	1000-1499		Set-Up	Units/Case Wgt.	sm@rt code	
	Includes: Post-it® DS330 Pop-Up Note Desktop Dispenser. Three (3) 3" x 3" 100-sheet Post-it® Pop-Up Note Pads in Canary Yellow with ruled lines packaged together in a Silver and Black Executive Gift Box.									6 sets/ 6.59 lbs.		
	No imprint	12.69	12.69	11.64	11.64	10.59	10.59	C	FREE			
	1 Color Ink Imprint on Dispenser Top (one or both sides of note opening)	N/A	18.18	15.62	13.10	11.66	11.43	CCCDD	FREE			
	2 Color Ink Imprint on Dispenser Top (one or both sides of note opening)	N/A	quote	quote	16.66	14.51	14.13	CDD	FREE			
	1 Laser Engraving* (one side of note opening)	N/A	17.88	16.46	15.71	13.47	13.34	CCCDD	FREE			
	Paper Insert, 1-4 color digital imprint**	N/A	19.48	15.84	14.15	12.91	12.21	CCCDD	FREE			
Supplied Business/Greeting Card Insert***	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE	FREE			

Allow 8 days production time. Actual Post-it® Note size: 2-7/8" x 2-7/8". Unprinted, 3" x 3" x 50- or 100-sheet Post-it® Pop-Up Note Pad refills are available retail. *Laser engraving is on underside of dispenser top on one side of opening. For laser engraving in more than one location, call for quote. For both laser engraving and ink imprint on same dispenser, call for quote and production time. **Color will shift slightly due to tinted top. Paper insert is digitally printed. All PMS inks will be converted for digital printing. Exact match to PMS ink color cannot be guaranteed. No fluorescent, special mix or metallic inks. *****Supplied insert cards must be sent within three (3) days of receipt of order.**

See pages 49-51 for product colors and standard ink imprint colors.

Conditions of sale page 46.



Custom Printed Products

Valuable education tools help create informative communications.


Education that Sticks

- Sheets
- Flags
- Page Markers
- Journals
- Posters

Post-it® Physician & Patient Education Tools allow you to provide medical professionals and patients with valuable educational materials.

With these customizable tools, you can help patients not only understand their disease and treatment, but also track their symptoms and navigate lifestyle changes. Physicians can benefit from a variety of drug information or exposure to new research data. Best of all, every product can stick practically anywhere.

Patient Education Ideas



Sheets can alternate with up to 4 different messages.

Alternating Designs, see page 37.

Take advantage of the ability to print cover sheets and up to four alternating inside sheets to make flexible **health tracking tools** that help patients monitor their disease to help maintain a healthy lifestyle.


Patient Compliance Ideas



Alternating Designs, see page 37.

Foster patients' **treatment compliance** with education. Using the cover as a training device, the back to inform and the inside sheet as a colorful take-home reminder.

Physician Education Ideas



Print each sheet with a different message.

Changeable Copy, see page 37.

Add a different one-color message on each sheet, creating a note tool that provides everything from **drug indications and contraindications** to dosage and side effects.


For more information on these versatile Post-it® Education Tools, please visit www.3M.com/ndl



Post-it® Flags on a Card, call for details.


Put information **where patients need it** with adhesive-backed posters that stick practically anywhere.

11" x 17"



Post-it® Note Pad attached to Poster Paper, call for details.


Multiple imprinting and customized design to create **appointment reminders** such as these, which help patients avoid scheduling conflicts by flagging upcoming appointments.




Custom-printed Kit, call for details.

Kits or Journals are available with any sheet or flag layout.

Track patient progress by providing a custom-printed kit, including handy forms that stick practically anywhere and flags to mark their next appointment.



Post-it® Poster Paper, see page 16.



Alternating Designs, see page 37.

Highlight **current research data** by using alternating sheets and four-color printing to create visually engaging notes that help complex statistics and data stick.

Announce **meetings and classes** with a temporary sign that sticks practically anywhere.



Custom Printed Products

We Offer Quick Response on Large Quantities!

For more information on these large quantity opportunities, please visit www.3M.com/ndl!



Post-it® Custom Printed Note and Flag Covers

Min. order: 10,000

Design your cover.

Combine features to fit your cover:

- 1" flag dispensers (TFD2D1)
- 1/2" duo flag dispensers (TFD2D2)
- 3" x 4" x 50-sheets
- 2" x 3" x 25-sheets
- 1" x 3" x 15-sheets
- Other sizes and options available, please call for details.

Custom print your cover and features in 1-4 colors!

Scotch®

NEW

Lint Sheet Pocket Packs — with 4-color Custom Printing

NEW full custom printing on large quantities! Custom print inside and outside of cover in 1-4 colors!



Post-it® Custom Printed Poster Paper

Vivid, photographic 1-4 color imagery on high-quality card stock paper. Full coverage Super Sticky adhesive makes mounting clean, crisp and neat. Choose 11" x 17" or 8-1/2" x 11" size. See page 16 for more product information.



NEW



Post-it® Custom Printed Notes Cubes in a Box

Min. order: 10,000

4 standard sizes available. Box is refillable. Custom print box in 1-4 colors with or without sheet printing. Large advertising area. Print a different message on each side. Custom dies available.

NEW

Post-it® Custom Printed Full Cover Gift Packs

Min. order: 10,000

NEW for Post-it® Flag Pens and Flag Highlighters. Double-pack holds 2 units, multi-pack holds 4 units. Custom print in 1-4 colors! Custom imprint barrels in up to 4 colors. Custom colors for barrels and flags also available.



Post-it® Custom Printed Designer 2 Dispensers

Min. order: 25,000

- Full customization!
- Imprint up to 4 colors inside dispenser on back card.
- Imprint up to 4 colors outside dispenser on top.
- Imprint flags in up to 2 colors (TFD2D1).



Post-it® Custom Printed Flag Card stock Dispensers

Min. order: 10,000

Choose 2-flag or 3-flag dispenser size. Custom print dispenser in 1-4 colors! Custom print flags in 1 or 2 colors!



Post-it® Custom Printed Flags in the Value Pak Dispenser

Min. order: 10,000

Up to 4 color imprint on the pack and 2 colors on the flags. Choose from a variety of flag colors. Low-profile package is ideal for mailing.




Post-it® Custom Printed DS330 Pop-Up Note Desktop Dispenser with Lenticular Insert Card

Min. order: 1,000

"WOW!" them with this eye-catching 3-D optical effect! Imprint lenticular insert card in up to 4 colors. Other imprint options available, for the DS330, see pages 28 and 29.

Great for "Before and After", Motion Graphics, Animations, Bold 3-D Impact Graphics!



Images animate back and forth as your customer holds the dispenser or walks up to use it!

Many Custom Options Available at Large Quantities

Embossing Foil Stamp Custom Sizes Custom Flag Highlighter Barrel and Flag Colors
Custom Shapes Custom Flag Pen Barrel and Flag Colors Custom Back Cards

Call for production time and quote.



Custom Printed
Notes

More Solutions to Fit Your Needs

Also Available...

For details on these products and more, please visit www.3M.com/ndl

Post-it® Custom Printed Notes Low Minimums

7 days production time

Call for production time on orders with more than 50 sets.

Available in three sizes: 3" x 4", 4" x 6" or 2" x 3-1/2" business card size. 50 sheets per pad. Your choice of standard White or Canary Yellow paper. 1-color imprint: choose standard Red, Blue or Black ink. Minimum order: 20 pads per design. Call for restrictions and details.



Custom Printed
Notes



7 days production time

Post-it® Super Sticky Notes are great for hard-to-stick surfaces. Notes are repositionable. Imprint 1-4 colors. White paper only. Available for all standard pad sizes. No imprint in adhesive area.

Visit www.3M.com/ndl for PDF brochure and details.



Post-it® Custom Printed Notes Executive Style Notes

7 days production time

Add 5 working days for covers.

Upscale for a professional impression. Choose from three elegant background designs: Granite, Linen or Parchment. Available in our most popular sizes. Custom imprint your logo or message. Add one additional color for background designs. Attractive protective cover available for 3" x 4" note size. Post-it® Flag Pens sold separately.



Create informative and effective communication and education tools with Post-it® Notes and Notes Cubes by using a variety of content options.

Great for promoting education programs, multiple services, brands, locations, products, mission statements, motion graphics, daily reminders and more!



Up to FOUR DIFFERENT 4-color designs alternate in one pad!

Alternating Designs

Post-it® Custom Printed Notes with **Alternating Designs**

Call for pricing and production time

Complete content change. . .print up to FOUR different 1–4 color designs to alternate through out entire pad. Popular note and notes cube sizes and sheet counts available.

Changeable Copy

Add **Changeable Copy** to your Post-it® Custom Printed Notes!

Call for pricing and production time

Ability to print a different 1 color (black) message on each sheet in addition to your 1–4 color note design. You choose the location of the message for each sheet!

Personalize for multiple customers in one order. Full-service mailable option available. Popular note and notes cube sizes and sheet counts available.



Add a different 1 color (black ink) message on each sheet of entire note pad!



Same 1–4 color design on every sheet.



Custom Printed Products

More Solutions to Fit Your Needs

Also Available. . .

For details on these products and more, please visit

www.3M.com/ndl

Post-it® Custom Printed Pop 'n Jot™ Refillable Small Cube Paks

8 days production time


Add 5 working days for sheet printing.

Imprint dispenser top and sides with 1–4 colors. Add impact with printed notes. Your choice of White or Canary Yellow paper.

Paper refills available.

**One-Hand Dispensing!
Stay-In-Place Base —
FREE UPON REQUEST!**

Notes automatically POP-UP one at a time!



Add your logo to the dispenser and flag pen!

Post-it® Custom Printed PRO330 Pop-Up Note Desktop Dispenser Gift Set

5 days production time

Stand up for the cause and personalize your Gift Set with your own logo or greeting. Each gift set includes: One (1) Post-it® PRO330 Pop-Up Note Desktop Dispenser, two (2) Pink 50-sheet Pop-Up Note Die-Cut Pads (flower shape) and one (1) Flag Pen with Pink Ribbon Imprint.

Plus send us your business or greeting card to be inserted at no extra charge.



Create Awareness

