

CLOSE OUT!

SPECIAL OFFER!

Valid Through December 31, 2009
or While Supplies Last!



Custom Printed
Flag Dispenser

*Organize and
Emphasize!*



Post-it® Mobile Flag Dispenser

- Loaded with 50 Post-it® Flags!
- Perfect for keychains, handbags, backpacks and anyone who's ON-THE-GO!
- Rotating cover protects the flags.

**From as low as
35¢ EACH!**

Post-it® Mobile Flag Dispenser		Quantity Price Per Unit						
Item #	Description and Imprint Options	100-249	250-499	500-999	1000-2499	2500-4999		Set-up
TFMD	Each Post-it® Mobile Flag Dispenser Includes: 50 Post-it® Flags, unprinted (Red only). White barrel with Silver metal keyring. Diameter: 5/8", Length: 3 5/8", Imprint Area: 1 1/2" x 3/8"							
	No Imprint (on barrel)	.35	.35	.35	.35	.35		C FREE
	1 color Imprint (on barrel)	Quote .94	1.16 1.51	.99 1.18	.89 1.04	.79 .92		C FREE
	2 color Imprint (on barrel) add	Quote .45	.10 .35	.08 .33	.06 .31	.04 .29		G FREE
Allow 5 days production time. Offer good only for products listed. Red flags only – flags cannot be customized. Individually wrapped in clear plastic sleeves and bulk packaged.				Be sure to reference offer #P3019 on purchase order				

**5 Day
Production
Time!**

Exclusive offers valid through December 31, 2009. While supplies last.