



Custom Printed Flag
Writing Instruments

NEW

Contour Series
Post-it® Flag
Highlighters, Pens & Markers

Each features contoured
comfort grip and integrated
Post-it® Flag dispenser!



highlight. write. flag. save.

Post-it® Flag Highlighters, Pens & Markers provide 2 tools in 1, so you can find what matters, fast!
2-in-1 convenience keeps your logo and message at your customer's fingertips!

Post-it® Flag Writing Instruments — Quantity Price Per Unit — Promotional Pricing Through December 31, 2009



NEW

Contour Series	NEW sleek design and comfort grip. 2-in-1 marking tool! Integrated Post-it® Flag dispenser comes loaded with 50 flags.	Item #	Imprint	100	250	500	1000	Set-Up	Be sure to reference offer #P4077 on your P.O.
Flag Pen (WIPFCCR) or Flag Highlighter (WIHFCCR) or Flag Permanent Marker (WIMFC)	1 Spot Color on Barrel			3.29	2.99	2.75	2.25	BBBC FREE	
	2 Spot Color on Barrel			4.39	3.49	3.25	2.55	BBBC FREE	
	3 & 4 Spot Color on Barrel			quote	4.19	3.75	2.85	BBC FREE	

Allow 5 days production time. No additional charge for PMS inks except fluorescent, special mixes and metallic inks \$285.00 (G). 4 color pricing does not include 4-color process printing. Flags cannot be customized. Pigment/ink is not refillable. Flag refills available with accompanying order, call for quote. Lanyards not available. Diameter: 9/16"; Length: 5-11/16"; Imprint Area on Barrel: 1-1/4" x 3/8"; Flags: 5/16", unprinted

Post-it® Flag Writing Instruments — Quantity Price Per Unit — Promotional Pricing Through December 31, 2009



FREE Lanyards
While Supplies Last!*

Industry Award Winning Product!

Trio Series	Post-it® Flag, Pen and Highlighter Combos provide 3 tools in 1 — highlight, write, and flag! Integrated Post-it® Flag dispenser comes loaded with 50 flags.	Item #	Imprint	100	250	500	1000	Set-Up	Be sure to reference offer #P4078 on your P.O.
Flag, Pen and Highlighter Combo (WIPHFC)	1 Spot Color on Barrel			3.58	2.83	2.59	2.24	BBBC FREE	
	2 Spot Color on Barrel			4.49	3.46	3.05	2.60	BBBC FREE	
	3 & 4 Spot Color on Barrel			quote	4.28	3.58	2.96	BBC FREE	
WILANYARDC: Lanyard — White Cord*				FREE	FREE	FREE	FREE	FREE	

Classic Series Post-it® Flag Highlighters, Pen & Markers provide 2 tools in 1, so you can find what matters, fast! Integrated Post-it® Flag dispenser comes loaded with 50 flags.

Classic Series	Post-it® Flag Highlighters, Pen & Markers provide 2 tools in 1, so you can find what matters, fast! Integrated Post-it® Flag dispenser comes loaded with 50 flags.	Item #	Imprint	100	250	500	1000	Set-Up	Be sure to reference offer #P4079 on your P.O.
Flag Pen (WIPFC) or Flag Highlighter (WIHF)	1 Spot Color on Barrel			3.29	2.50	2.20	1.96	BBBC FREE	
	2 Spot Color on Barrel			4.10	3.12	2.72	2.25	BBBC FREE	
	3 & 4 Spot Color on Barrel			quote	3.89	3.19	2.60	BBC FREE	
WILANYARD: Lanyard — Choose White, Yellow, Blue or Black Cord*				FREE	FREE	FREE	FREE	FREE	

Allow 5 days production time. No additional charge for PMS inks except fluorescent, special mixes and metallic inks \$285.00 (G). 4 color pricing does not include 4-color process printing. Flags cannot be customized. Pigment/ink is not refillable. Flag refills available with accompanying order, call for quote.

*FREE Lanyards While Supplies Last. Lanyards available upon request only with accompanying order. Please indicate color choice for WIPFC or WIHF (White, Yellow, Blue or Black cord) at time of order. White lanyards only for WIPHFC. Not available for Contour Series.



Special Offers Valid September 1 – December 31, 2009

3M and Post-it are trademarks of 3M. © 3M 2009. All Rights Reserved.

