



Custom Printed
Products

Health Care

Health care professionals rely on 3M and Post-it® Products every day. Get your product, brand or message repeated exposure with custom printed Post-it® Notes, Post-it® Flag Pens and Highlighters or a quality 3M™ Littmann® Stethoscope.

Post-it® Die Cut Notes

Post-it® Die Cut Notes help you stand out and communicate your message with shapes that enhance your brand and get their attention.

Post-it® Super Sticky Notes

Post-it® Super Sticky Notes hold stronger and longer than most sticky notes and are great for vertical surfaces like medical equipment and exam room doors so your message gets noticed again and again.

Post-it® Flag Highlighters and Pens with Lanyard

Put your message on the product that lets health care workers highlight or notate a section, then flag the page for later reference. They're great for working with medical charts and records.

Post-it® Note Cube

Long lasting impressions with side-printing that is easily visible on the desktop and individually printed notes that carry your message throughout the office with each sheet.

3M™ Littmann® Stethoscope

A quality 3M™ Littmann® Stethoscope bearing your logo is a powerful way to make a lasting and continuing impression on your physician audience.

Who
do
you
want
to
reach?



The brand they trust,
the brand you rely on.

3M, Littmann, Post-it and the color Canary Yellow are trademarks of 3M.

© 3M 2007