

## Exact Color Match



Custom Printed  
Products

*Gift Sets*

### 5–8 days production time

- Available undecorated or with a custom imprint on the Post-it® Product.
- Gift set comes packaged in an Executive Gift Box.
- Send us your business card or personalized greeting card to be inserted at no extra charge!† See page 27 for availability.
- See page 43 for available colors.

A professional gift to show your appreciation—anytime of the year!



Dispenser comes boxed with three Post-it® Pop-Up Note Pads!

TPJDS330GB



TPJDS100GB




Mix and match any two Classic Series Post-it® Flag Writing Tools!



WIPHGB

**Post-it® Products Executive Gift Sets**

**Quantity Price per Gift Set**

	Item #	Imprint	100-	250-	500-	1000-	2500-	5000-		Set-up	Sets/ Case Wt.	3M@r	
			249	499	999	2499	4999	9999					
<b>Writing Tool Executive Gift Set</b> 	WIPHGB	1 Spot Color Imprint (On Barrel of Both Items)	7.20	6.05	5.38	4.90	4.60	4.05	BBBBCCD	FREE	250/ 25 lbs.	4S	
		2 Spot Color Imprint (On Barrel of Both Items)	8.30	7.03	6.26	5.66	5.30	4.27	BBBBCCD	FREE			
		3 & 4 Spot Color Imprint (On Barrel of Both Items)	quote	quote	7.46	6.76	6.03	4.59	BCCD	FREE			
		<p><b>Allow 5 days production time.</b>                      Includes: Any one (1) Classic Series Post-it® Flag Pen WIPFC and any one (1) Flag Highlighter WIHFC packaged together in a Silver Gift Box. Each writing tool holds 50 flags. Imprint your logo or message on barrel of both items — same art on both.</p> <p>Writing tool ink and pigment are not refillable. Flags cannot be customized. 4 color pricing does not include 4-color process printing. Trio Series and Contour Series not available for this gift set.</p>											
<b>DS330 Pop-Up Note Dispenser Executive Gift Set*</b> 	TPJDS330GB	No Imprint	12.69	12.69	11.64	11.64	10.59	10.59	C	FREE	6/ 6.59 lbs.	4S 4S 1X 4C	
		1 Spot Color Imprint on Dispenser Top (One or both sides of opening)		18.18	15.62	13.10	11.66	11.43	CCCCD	FREE			
		2 Spot Color Imprint on Dispenser Top (One or both sides of opening)	quote	quote	16.66	14.51	14.13	CDD	FREE				
		1 Laser Engraving on Dispenser Top** (One side of opening)		17.88	16.46	15.71	13.47	13.34	CCCCD	FREE			
		Paper Insert, 1-4 Color Digital Print***		19.48	15.84	14.15	12.91	12.21	CCCCD	FREE			
		Supplied Business/Greeting Card Insert†	FREE	FREE	FREE	FREE	FREE	FREE					
		<p><b>Allow 8 days production time.</b>                      Includes: One (1) DS330 Pop-Up Note Dispenser and three (3) 3" x 3" 100-sheet Post-it® Pop-Up Note Pads* in Canary Yellow with ruled lines packaged together in Silver and Black Executive Gift Box.</p>											
<b>DS100 Pop-Up Note Dispenser Executive Gift Set*</b> 	TPJDS100GB	No Imprint	17.94	17.94	17.94	16.89	16.89	15.84	C	FREE	10/ 17.44 lbs.	4S 4S 1X	
		1 Spot Color Imprint on Dispenser Top (One or both sides of opening)			23.08	20.86	18.46	16.88	16.67	CCCCD			FREE
		2 Spot Color Imprint on Dispenser Top (One or both sides of opening)			26.76	24.01	20.98	18.87	18.56	CCCCD			FREE
		1 Laser Engraving on Dispenser Top** (One side of opening)			23.43	21.93	20.13	18.56	16.96	CCCCD			FREE
		Supplied Business/Greeting Card Insert†	FREE	FREE	FREE	FREE	FREE	FREE	FREE				
		<p><b>Allow 8 days production time.</b>                      Includes: One (1) DS100 Pop-Up Note and Flag Dispenser; one (1) 3" x 3" 50-sheet Post-it® Pop-Up Note Pad* in Canary Yellow (unprinted); one (1) Red 1" Post-it® Flag Dispenser and one (1) Blue 1" Post-it® Flag Dispenser (100 flags total) packaged together in Silver and Black Executive Gift Box.</p>											

\*Actual note size is 2-7/8" x 2-7/8". Unprinted, 3" x 3" 50- or 100-sheet Post-it® Pop-Up Note Pads and Post-it® Flag Dispensers are available retail.

\*\*Laser engraving is on underside of dispenser top on one side of note opening. For both laser engraving and ink imprint on the same dispenser or for laser engraving in more than one location, call for quote and production time.

\*\*\*Paper insert is digitally printed. All PMS inks will be converted to digital printing. Exact match to PMS ink color cannot be guaranteed.

†Supplied insert cards must be sent within three (3) days of receipt of order.

Conditions of Sale page 38.

Gift Sets

## Create Awareness!



Stand up for the cause and personalize your Gift Set with your logo!

**Post-it® Custom Printed PRO330 Pop-Up Note Dispenser Gift Set**

Each gift set includes: One (1) Post-it® PRO330 Pop-Up Note Dispenser, two (2) Pink 50-sheet Pop-Up Note Die-Cut Pads (flower shape) and one (1) Flag Pen with Pink Ribbon Imprint. Plus send us your business or greeting to be inserted† at no extra charge.

**5 Days Production.**  
 For more details, visit:  
[www.3M.com/ndl](http://www.3M.com/ndl)



Custom Printed Products

Call for production time and pricing.

**We Offer Quick Response on Large Quantities!**

*Many Custom Options Available at Large Quantities:*

- Embossing
- Foil Stamp
- Custom Sizes
- Custom Barrel Colors
- Custom Shapes
- Custom Flag Colors
- Custom Back Cards

For more information on these large quantity opportunities please visit [www.3M.com/ndl](http://www.3M.com/ndl)



**Post-it® Custom Printed Note and Flag Covers**

**Min. order: 10,000**

Design your cover. Combine features to fit your cover:

- 1" flag dispensers (TFD2D1)
- 1/2" duo flag dispensers (TFD2D2)
- 3" x 4" x 50-sheets
- 2" x 3" x 25-sheets
- 1" x 3" x 15-sheets
- Other sizes and options available.



**Post-it® Custom Printed Notes Cube in a Box**

**Min. order: 10,000**

4 standard sizes available. Box is refillable. Custom imprint box in up to 4 colors with or without sheet printing. Large advertising area. Print a different message on each side. Custom die-cut backers available.



**Post-it® Custom Printed Full Cover Gift Packs**

**Min. order: 10,000**

Classic Series Post-it® Flag Writing Tools. Double-pack holds 2 tools, multi-pack holds 4 tools. Custom imprint cover in up to 4 colors! Custom imprint barrels in up to 4 colors. Custom colors for barrels and flags also available.



**Post-it® Custom Printed Flag Card Stock Dispensers**

**Min. order: 10,000**

Choose 2-flag or 3-flag dispenser size. Custom imprint dispenser in up to 4 colors! Custom imprint flags in 1 or 2 colors!



**Post-it® Custom Printed Designer 2 Dispensers**

**Min. order: 25,000**

Full customization — inside and out!

- Imprint up to 4 colors inside dispenser on back card.
- Imprint up to 4 colors outside dispenser on top.
- Imprint flags in 1 or 2 colors (TFD2D1 only).



**Post-it® Custom Printed Flags in the Value Pak Dispenser**

**Min. order: 10,000**

Custom imprint up to 4 colors on the pak and 2 colors on the flags. Choose from a variety of flag colors. Low-profile package is ideal for mailing.

# More Solutions to Fit Your Needs

★ Larger Quantity Opportunities More Solutions...

Create informative and effective communication and education tools with Post-it® Notes and Notes Cubes by using a variety of content options.

Great for promoting education programs, multiple services, brands, locations, products, mission statements, motion graphics, daily reminders and more!

Up to FOUR DIFFERENT 4 color designs alternate in one note pad!

## Post-it® Custom Printed Notes with Alternating Designs

### Call for pricing and production time

Complete content change...custom print up to FOUR DIFFERENT 1 to 4 color designs to alternate through out entire note pad. Popular note and notes cube sizes and sheet counts available.

Sheets can alternate with up to 4 different designs.

Add a different 1 color (black ink) message on each sheet of entire note pad. PLUS same 1 to 4 color design on every sheet!

## Add Changeable Copy to your Post-it® Custom Printed Notes!

### Call for pricing and production time

Ability to custom print a different 1 color (black) message on each sheet in addition to your 1 to 4 color note design. You choose the location of the message for each sheet!

Personalize for multiple customers in one order. Full-service mailable option available. Popular note and notes cube sizes and sheet counts available.

Print each sheet with a different message.



Custom Printed Products

**Call for restrictions and pricing.**

*Customers not only love Post-it® Notes, but are conditioned to use them over and over again... which keeps your message visible day after day and sheet after sheet!*

For more information on these products and more please visit [www.3M.com/ndl](http://www.3M.com/ndl)



**Quantities under 500!**

**Post-it® Custom Printed Notes Low Minimums**

**7 days production time**

**Call for production time on orders with more than 50 sets.**

Available in three sizes: 3" x 4", 4" x 6" or 2" x 3-1/2" business card size. 50 sheets per note pad. Your choice of standard paper: White or Canary Yellow. 1 color imprint: choose standard Red, Blue or Black ink. Minimum order: 20 note pads per design.



Custom Printed Notes




**Post-it® Custom Printed Super Sticky Notes**

**7 days production time**

Post-it® Super Sticky Notes are great for hard-to-stick surfaces. Notes are repositionable. Custom imprint in up to 4 colors. No imprint in adhesive area.

**Available for all standard note sizes! White Paper Only! 10% Upcharge. Visit [www.3M.com/ndl](http://www.3M.com/ndl) for details and restrictions.**



**Designed for vertical surfaces!**



**Post-it® Custom Printed Flag Dispenser**

**5-8 days production time**

Large area for up to 4 color custom imprint.\* Choose from a variety of flag colors. Holds 50 1" flags.

\*PMS colors will be converted to the closest possible 4-color process match. Exact color match cannot be guaranteed. No fluorescent, special mix or metallic inks.

**Post-it® Custom Printed Backsheet**

**Add 1 day to normal note pad production time**

*Backsheets = liner of note pad.*  
 Backsheets can be printed as a business reply card, customer survey, coupon or whatever you choose. Makes a great final reminder to complement your other message. Popular sizes available.



4" x 6" size shown

**Post-it® Custom Printed Pop 'n Jot™ Refillable Small Cube Packs**

**8 days production time**

**Add 5 working days for sheet printing.**

Imprint dispenser top and sides with up to 4 colors. Add impact with printed notes. Your choice of White or Canary Yellow paper. Paper refills available.



**One-Hand Dispensing!  
 Stay-In-Place Base —  
 FREE UPON REQUEST!**

Notes automatically POP-UP one at a time!



**3M™ Littmann® Stethoscopes**

**Call for production time**

The brand of choice for over 40 years! Engrave your message on the side of the binaurals, the yoke area or tubing. Available in three models and a variety of tubing colors.



**Customer Solutions**



Trade Shows



Medical/Healthcare



Finance



Business Professionals



Education



Forms & Applications