

← Minimum Bleed Margins  
 ← Final Pad Size  
 ← Imprint Margins  
 ← Fold Lines

For custom printed notes, see next page.

- RESTRICTIONS**
- Logos, text or anything that needs to be guaranteed to print should be placed within the imprint margin.
  - Bleeds should fill the minimum bleed margin.
  - Convert all imprint colors to CMYK (no spot or RGB colors).
  - Bleeds allowed on minimum quantity orders of 250.

Approximate Final size: 4.125"W x 8.218"H

Bleed Area: 4.375"W x 8.4687"H


Outside Cover Imprint Areas (denoted by black dotted lines):

Area A: 3.875"W x 1.757"H

Area B: 3.875"W x 2.843"H

Area C: 3.875"W x 0.642"H

Area D: 3.875"W x 0.6"H



**RESTRICTIONS**

- Keep text at/or above 4 pt. and not too thin.
- TM and ® — minimum of 3 pt.
- Line thickness — minimum of .3 pt.
- Bleeds should fill the minimum bleed margin.
- Logos, text or anything that needs to be guaranteed to print should be placed within the imprint margin.
- Convert all imprint colors to CMYK (no spot or RGB colors).

Approximate Final Pad Size: 4"W x 2.875"H

Approximate Imprint Area: 3.625"W x 2.5"H

Approximate Bleed Area: 4.25"W x 3.125"H

Beneath each sheet is adhesive area